

OSYPKA iSense[®]

ELECTRICAL CARDIOMETRY™ SENSORS




Window to the Circulation®

OSYPKA iSense® sensors are intended for use with the Electrical Cardiometry™ (EC™) monitor AESCULON® and cardiometers ICON® and ICONCore®. OSYPKA iSense® Sensors are 100% disposable and are available in various sizes, for neonate, pediatric, and adults. All sensors utilize a highly conductive AgCl wet gel and strong hypoallergenic adhesion, which provides instant electrical contact for high quality signal conductivity. Four iSense® sensors are packed together in one pouch for the convenience of the user and more importantly keeping the gel fresh for best conductivity and signal quality, specially in difficult patients.



OSYPKA iSense® sensor box



Patient cable

OSYPKA iSense® Specification:

Sensor	Patient Group	Type/ Size	Article No.	Packaging Unit
	Premature and Pre-mature	Prewired sensor 22 mm x 22 mm (rectangular)	XN.TSP	20 pouches
	Adults and Children	Sensor with clip 30 mm (round)	XN.TSN	25 pouches

Products and options may not be available in every region or country. Please contact mail@osypkamed.com for availability in your region and country and for further information.

Osypka Medical GmbH

Albert-Einstein-Strasse 3
12489 Berlin, Germany
Phone: +49 (30) 6392 8300
Fax: +49 (30) 6392 8301
E-Mail: mail@osypkamed.com
www.osypkamed.com

United States of America: Osypka Medical, Inc.

dba. Cardiotronic
La Jolla, CA 92037, USA
Phone: +1 (858) 454 0600
Fax: +1 (858) 454 0640
E-Mail: mail@cardiotronic.net
www.cardiotronic.net

U.S. Patent Nr. 6,511,438. Other patents pending.

iSense, AESCULON, ICON, ICONCore, Cardiotronic, Electrical Cardiometry, Window to the Circulation and logos are trademarks of Osypka Medical.

Copyright © 2024 Osypka Medical

